

Diagnostic Display



| Product Model | C41W+ | C42W+ | C61W+ | C81W+ |
|------------------------|--|--|---|---|
| Product Name | 4MP Integrated Dual-screen Display | 4MP Integrated Dual-screen Display | 6MP Integrated Dual-screen Display | 8MP Integrated Dual-screen Display |
| Display size | 27" | 30" | 30" | 31.1" |
| Support color | 1.07 Billion | 1.07 Billion | 1.07 Billion | 1.07 Billion |
| Resolution | 2560 x 1440 | 2560 x 1600 | 3280x2048 | 4096 x 2160 |
| View angle | 178° Typ. | 178° Typ. | 170° Typ. | 176° Typ. |
| Brightness | 350cd/m² Typ. | 700cd/m² Typ. | 800cd/m² Typ. | 850cd/m² Typ. |
| Contrast | 1000: 1 Typ. | 1000: 1 Typ. | 1000: 1 Typ. | 1450: 1 Typ. |
| Input signal interface | DVI-D : 24 pin * 2 DisplayPort : 20 pin DIP:19 pin | DVI-D : 24 pin * 2 DisplayPort : 20 pin DIP:19 pin | DVI- D: 24 pin * 2 DisplayPort: 20 pin * 2 | DVI- D: 24 pin * 2 DisplayPort: 20 pin * 2 |



| Product Model | C22S+ | C22SP+ | C24S+ | G22S+ | G22SP+ | C32S+ | C32SP+ | G32S+ | G32SP+ | G52S+ | G52SP+ | |
|------------------------|---|---|---|---|---|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|
| Product Name | 2MP Color Medical Display | 2MP Color Medical Display | 2MP Color Medical Display | 2MP Grayscale Medical Display | 2MP Grayscale Medical Display | 3MP Color Medical Display | 3MP Color Medical Display | 3MP Grayscale Medical Display | 5MP Grayscale Medical Display | 5MP Grayscale Medical Display | 5MP Grayscale Medical Display | |
| Display size | 21.3" | 21.3" | 21.3" | 21.3" | 21.3" | 21.3" | 21.3" | 21.3" | 21.3" | 21.3" | 21.3" | |
| Support color | 1.07 Billion | 1.07 Billion | 1.07 Billion | 4096 Grayscale | 4096 Grayscale | 1.07 Billion | 1.07 Billion | 4096 Grayscale | 4096 Grayscale | 4096 Grayscale | 4096 Grayscale | |
| Resolution | 1600 x 1200 | 1600 x 1200 | 1600 x 1200 | 1600 x 1200 | 1600 x 1200 | 2048 x 1536 | 2048 x 1536 | 2048 x 1536 | 2048 x 1536 | 2560 x 2048 | 2560 x 2048 | |
| View angle | 176° Typ. | 176° Typ. | 178° Typ. | 176° Typ. | 176° Typ. | 176° Typ. | 176° Typ. | 176° Typ. | 176° Typ. | 170° Typ. | 170° Typ. | |
| Brightness | 900cd/m² Typ. | 900cd/m² Typ. | 440cd/m² Typ. | 1900cd/m² Typ. | 1900cd/m² Typ. | 800cd/m² Typ. | 800cd/m² Typ. | 1700cd/m² Typ. | 1700cd/m² Typ. | 1200cd/m² Typ. | 1200cd/m² Typ. | |
| Contrast | 1400: 1 Typ. | 1400: 1 Typ. | 1500: 1 Typ. | 1400: 1 Typ. | 1400: 1 Typ. | 1400: 1 Typ. | 1400: 1 Typ. | 1400: 1 Typ. | 1400: 1 Typ. | 1200: 1 Typ. | 1200: 1 Typ. | |
| Input signal interface | DVI- D: 24 pin VGA: D- Sub 15 pin DisplayPort: 20 pin | DVI- D: 24 pin VGA: D- Sub 15 pin DisplayPort: 20 pin | DVI- D: 24 pin VGA: D- Sub 15 pin DisplayPort: 20 pin | DVI- D: 24 pin VGA: D- Sub 15 pin DisplayPort: 20 pin | DVI- D: 24 pin VGA: D- Sub 15 pin DisplayPort: 20 pin | DVI- D: 24 pin DisplayPort: 20 pin | DVI- D: 24 pin DisplayPort: 20 pin | DVI- D: 24 pin DisplayPort: 20 pin | DVI- D: 24 pin DisplayPort: 20 pin | DVI- D: 24 pin DisplayPort: 20 pin | DVI- D: 24 pin DisplayPort: 20 pin | DVI- D: 24 pin DisplayPort: 20 pin |