

# GIIS

## 1MP Grayscale Medical Display



## 1MP Grayscale Medical Display

With built-in multiple display mode, GIIS can satisfy doctors' requirement of reviewing multiple images. With excellent display performance, GIIS provide the perfect clinical image for hospital, help doctors to make clinical decisions quickly, improving hospital efficiency and patient satisfaction.

Each monitor is strictly calibrated to compliant with DICOM standard, ensure the accurate display of medical images and details. The built-in backlight stabilization system can ensure the consistence performance in long-term use.

# BEACON

## 1MP Grayscale Medical Display

### Application Area

- DSA
- DSI
- CT
- MRI



### Technical Specification



Built-in multiple display mode



Backlight stabilization system



Compliant with DICOM standard



Wide viewing angel

Model	G11S
Display size	19"
Support color	1024 grayscale
Resolution	1280 x 1024
Pixel Pitch	0.294 x 0.294 mm
Brightness	1300 cd/m2 Typ.
Contrast	1000:1 Typ.
Angle of view	178° Typ.
Response time	25 ms Typ.
Input signal interface	DVI: DVI-D 24pin VGA: D-sub 15pin Video(SOG): BNC IN BNC OUT
Size	425(W)x449~523(H)x225(D) mm
Weight	6.3kg
Standard Mounting	VESA 100x100mm
Power supply	AC 100~240V
Power Consumption	55W max.
Power Saving	EPA, VESA
MTBF	> 50,000 Hours