

# CI4S

## 1MP Color Medical Display



## 1MP Color Medical Display

With built-in multiple display mode, CI4S can satisfy doctors' requirement of reviewing multiple images. With excellent display performance, CI4S provide the perfect clinical image for hospital, help doctors to make clinical decisions quickly, improving hospital efficiency and patient satisfaction.

Each monitor is strictly calibrated to compliant with DICOM standard, ensure the accurate display of medical images and details. The built-in backlight stabilization system can ensure the consistence performance in long-term use.

# BEACON

## 1MP Color Medical Display

### Application Area

- DSA
- DSI
- CT
- MRI



### Technical Specification



Built-in multiple display mode



Backlight stabilization system



Compliant with DICOM standard



Wide viewing angel

Model	CI4S
Display size	19"
Support color	16.7M
Resolution	1280 x 1024
Pixel Pitch	0.294 x 0.294 mm
Brightness	300 cd/m2 Typ.
Contrast	1000:1 Typ.
Angle of view	178° Typ.
Response time	14 ms Typ.
Input signal interface	DisplayPort: 20pin DVI: DVI-D 24pin VGA: D-sub 15pin
Size	425(W)x449~523(H)x225(D) mm
Weight	6.3kg
Standard Mounting	VESA 100x100mm
Power supply	AC 100~240V
Power Consumption	50W max.
Power Saving	EPA, VESA
MTBF	> 50,000 Hours