


G32S+

3MP Grayscale Medical Display



3MP Grayscale Medical Display

G32S+ is with 3MP (2048X1536) high resolution, high brightness, high contrast ratio and wide viewing angle, can be applied to various medical image equipment and PACS system.

Each monitor is strictly calibrated to compliant with DICOM standard, ensure the accurate display of medical images and details. The built-in backlight stabilization system can ensure the consistence performance in long-term use.

G32S+ is with LED backlight and ultra-thin design.

BEACON

Application Area

- CR/DR
- CT
- MRI
- X-RAY
- DSA
- DSI



Technical Specification



Built-in multiple display mode



Backlight stabilization system



Ambient light adaptive system(ALAS)



QA software



Compliant with DICOM standard

Model	G32S +
Display size	21.3 "
Support color	4096 grayscale
Resolution	2048x1536
Pixel Pitch	0.2115x0.2115mm
Brightness	1700cd/m2 Typ.
Contrast	1400:1 Typ.
Angle of view	176° Typ.
Response time	20ms Typ.
Input signal interface	DVI-D: 24 pin Display Port:20 pin
Size	369(W)x507~627(H)x235(D)mm
Weight	9.2 kg
Standard Mounting	VESA 100x100mm
Power supply	DC 12V
Power Consumption	37W max.
Power Saving	EPA, VESA
MTBF	> 50,000 hours