

C22W

21.5" Multi-function Clinical Display



21.5" Multi-function Clinical Display

Fully compliant with the DICOM standard and with FHD widescreen design, C22W can help doctors to access to patient information in a filmless environment while browsing different medical images. That can be easily operated, and improve doctors' working efficiency and patient satisfaction.


BEACON

C22W

21.5" Multi-function Clinical Display

Application Area

- Outpatient Service
- Sickroom
- Intensive care
- Clinical department
- Community medical system



Technical Specification



Compliant with DICOM standard



Built-in multiple display mode



Wide view angle

Model	C22W
Display size	21.5 "
Support color	16.7M
Resolution	1920x1080
Pixel Pitch	0.248x0.248mm
Brightness	250 cd/m2 Typ
Contrast	1000:1 Typ
Viewing angle	Horizontal 178°, vertical 178° (CR > 10)
Response time	14ms Typ .
Input signal interface	DisplayPort: 20pin
	DVI: DVI-D24pin
	VGA: D-sub 15pin
	Audio: 3.5mm stereo jack
Size	Speaker: 5W*2
	530(W) x 320(H) x 70 (D) mm
	Weight
	4.5kg
Standard mounting	VESA 100x100mm
Power supply	AC 110-240 50/60HZ 1A
Power consumption	40W max
Power saving	EPA,VESA
MTBF	> 30,000 Hours

BEACON

Shenzhen Beacon Display Technology Co.,Ltd.

12F, Block B1, Nanshan Zhiyuan, No.1001
Xueyuan Road, Nanshan District,
Shenzhen, 518055 China

Tel: +86 755 2699 5355
Fax: +86 755 2699 5755
www.beacon-display.cn

All specification are subjected to change without prior notice.