



## PRODUCT FEATURES

### SPIROMETRY FILTER 59 "OVAL"

*For use during pulmonary lung function testing*

**Description:** Barrier filtration for protection of patient, staff and equipment from bacterial and viral contamination.  
Filter with "oval semi-rigid anatomical mouthpiece" to allow adequate blockage of mouth while avoiding possible leaks and dispersion from the sides of the mouth itself.

**Code items:**

- individually packaged:

**30803 ADULT**  
**30804 PAEDIATRIC**

**Materials:**

- housing:  
- filtration membrane:  
- mouthpiece

Polystyrene (HIPS) recyclable  
Hydrophobic Polypropylene (SMS)  
Pvc semi-rigid (non-toxic)

**Kind of filtration:**

Electrostatic / Mechanical

**Flow resistance: (complies with ATS)**

- flow rate 12,00 l/s 0,93 cm H<sub>2</sub>O / (l/s)

**Cross Contamination Reduction Efficiency:** > 99,99 %

**Dead space:** 40,53 cm<sup>3</sup>  
(connector not included)

**Filtration area:** 25,62 cm<sup>2</sup>  
**Exchange filtration area:** > 72,00 cm<sup>2</sup>

**Weight:** 22,00 g

**Working procedure:**

CE 0373 Class IIa - Italian Superior Institute of Health  
Disposable

**Quantity / box:** 75 pcs.

**NOTES:** single patient use. Latex free.

**Any size of adaptors is available upon request.**

**MADA**

*Spirometry Filters*