



Single Syringe Kits

Salient 190 mL Syringe and Quick Fill Tube



Overview


Catalogue number	ZY6322
Globan number	84710261
Quantity	50/box
Shelf life	5 years

Packaging

Primary	Tyvek 1073B & 0.066 mm PE/PA clear
Secondary	Cardboard shipper box
Assembled in	China

Product information

PSI	300 PSI
Volume	190 mL
Length	186.84 mm
Diameter internal	46.84 mm
Diameter external	50.85 mm; Flange: 64.01 mm
Sterilization method	Ethylene Oxide (EtO)



Component	Material
1 End Cap (MDJP015)	Huntsman P5M4K-046 Polypropylene
2 Plunger (MDAS001)	Makrolon 2530-451118
Rubber Seal (DJP002)	Santoprene TPV 181-55med
3 Syringe barrel	Eastman PET MN052
4 Dust Cap (MTJP015)	Huntsman P5M4K-046 Polypropylene
5 Quick Fill Tube (80105)	Dow #9931 Low density Polyethylene