

Product Summary Sheet

Catalog/Product/Reference Code: **150001**

Description: **T-type pressure tube with check valve**

Components/Material

Components		Material
Name	Qty	
Female luer Joint	2	Clear w/Polycarbonate(PC)
Male luer Joint	1	Clear PVC
Female Dust Cap	1	White HDPE
Male Dust Cap	2	Blue PP
Check valve	1	Medical Silicone Rubber
Tube	1	φ1.7*φ3.5*1500mm clear PVC

Packaging and Sterilization

Attribute	Detail
Inner package	60g/m ² medical use paper+ PET/LDPE film
Medium bag	White box
Outer carton	White carton

Technical Characteristics

Attribute	Detail
Pressure rating	350psi
Internal volume	3.4ml
Intended use	Special for double barrel and MRI syringe, To deliver, in low-pressure injection system, contrast media to the vascular system for diagnostic purposes

PHOTO

