

Product Summary Sheet

Catalog/Product/Reference Code: **150003**

Description: **T-type pressure tube with check valve**

Components/Material

Components		Material
Name	Qty	
Female luer Joint	2	Clear w/Polycarbonate(PC)
Male luer Joint	1	Clear PVC
Female Dust Cap	1	White HDPE
Male Dust Cap	2	Blue PP
Check valve	1	Medical Silicone Rubber
Tube	1	φ1.7*φ3.5*2500mm clear PVC

Packaging and Sterilization

Attribute	Detail
Inner package	60g/m ² medical use paper+ PET/LDPE film
Medium bag	White box
Outer carton	White carton

Technical Characteristics

Attribute	Detail
Pressure rating	350psi
Intended use	Special for double barrel and MRI syringe, To deliver, in low-pressure injection system, contrast media to the vascular system for diagnostic purposes

PHOTO



Label Representation

