

MRI Stethoscope



100% antimagnetic stethoscope with 56 cm long double tube (latex-free).

The two-tube system guides the sound waves separately to each ear, thus creating the conditions for perfect auscultation (listening to the body).

Double-head chest-piece with great advantage of bell mode for low-frequency sounds 20 Hz to 100 Hz and membrane for high frequencies 100 Hz to 1,000 Hz.

Double-head chest piece can also be rotated 180 degrees if required.

Three additional bells (adult, child and infant) so it can also be used as a child or infant stethoscope.

Three different ear tips for optimal fitting are included in the set.

Includes two extra membranes as replacements.

Excellent acoustic quality and versatility for daily use.