

# UP-D899MD

A6 Digital Black and White Thermal Printer



## Overview

The UP-D899MD is a compact A6 black and white digital printer and is the successor to Sony's medical printer UP-898 series. It delivers high speed, outstanding quality black & white printing with enhanced image accuracy, an improved user interface/printing cutter, and upgraded environmental performance.

## Features

### Wide PC compatibility

For maximum flexibility, a wide range of PC drivers are provided for Windows 10,11. When choosing compatible mode, you can select 898/897 driver.

### Compliance with medical standards

This product is distributed to the US and EU as a medical device and satisfies product safety standards (e.g., IEC 60601-1).  
For more details, please contact your nearest Sony sales office or an authorized dealer.

## Specifications

| Specifications       |  |
|----------------------|--|
| Printing method      | Direct thermal printing  |
| Resolution           | 325 dpi  |
| Gradation            | 8 bits (256 levels) processing   |
| Picture elements     | Digital: 1,280 x 2,560 dots  |
| Picture area         | 100 x 200 mm (3 7/8 x 7 7/8 inches)  |
| Paper size           | Paper width of 110 mm (4 3/8 inches)   |
| Printing time        | High-speed mode: Approx. 2 seconds/image (at standard setting) (1,280 x 960 dots)<br>Normal speed mode: Approx. 3.8 seconds/image (at standard setting) (1,280 x 960 dots) |
| Interface            | Hi-Speed USB connector (USB 2.0): type C (for data receiving) (x1)   |
| Printer Driver SW    | Microsoft Windows 10, 11   |
| Picture requirements | AC 100 V to 240 V, 50/60 Hz  |
| Input current        | 0.75 A to 0.5 A  |

## Specifications

|                                   |  |
|-----------------------------------|--|
| Operating temperature             | 10°C to 40°C (50°F to 104°F)   |
| Operating humidity                | 20% to 80% (no condensation allowed)   |
| Storage and transport temperature | -20°C to +60°C (-4 °F to +140°F)   |
| Storage and transport humidity    | 20% to 90% (no condensation allowed)   |
| Dimensions (W*H*D)                | 154 × 88 × 208 mm (6 1/8 × 3 1/2 × 8 1/4 inches)   |
| Mass                              | 2.0 kg (4 lb 6.5 oz) (printer only)  |
| Supplied accessories              | Thermal head cleaning sheet (1)<br>CD-ROM(1)<br>Before Using this Printer (1)<br>Service Contact List (1)<br>Information for Customers in Europe (1)   |
| EC directive                      | * The values for dimensions are approximate.<br>This product is distributed to the US and EU as a medical device and satisfies product safety standards (e.g. IEC 60601-1)**.<br>** For more details, please contact your nearest Sony sales office or an authorized dealer. |

Information on EU Data Act: [Click here](#) to check if your Sony product and/or related service is subject to the EU Data Act.

## Gallery





---

© 2004 - 2025 Sony Corporation. All rights reserved. Reproduction in whole or in part without written permission is prohibited. Features and specifications are subject to change without notice. The values for mass and dimension are approximate. All trademarks are the property of their respective owners.